



SURFIN BETA

REMOVER FOR ALBUMINOID STAINS: MILK, EGG, SWEAT, URINE, BLOOD, ICECREAM, VOMIT, WATER HALOES ETC.

SURFIN BETA is a combination of particular spotting substances with solvents and emulgators.

Aspect	Limpid liquid
Color	Transparent
Smell	Characteristic
Density at 20°C	1 g/ml
pH at 20°C	8
Freezing temperature	Less than -7°C
Solubility in water	Clear solution
Solubility in perchloroethylene	Dispersion

SURFIN BETA can be used both in drycleaning and laundering, for pre-spotting and post-spotting.

SURFIN BETA can be easily rinsed in perchloroethylene and water.

Pre-spotting:

Dampen the stain with SURFIN BETA and rub delicately off, buffer with a detergent (e.g. SURFAN) and then dryclean.

Post-spotting:

Dilute SURFIN BETA 1:5 with water, use water or steam spotting machines, and dry well with air.

For delicate fabrics and weak colours it is advisable to test the product on an internal sewing of the garment before use it.

SURFIN BETA can be stored for 3 years.

The bottle should always be tightly closed after the use.

Plastic bottle with 450 ml content.

Plastic container with 5 - 10 Kg content.

714250 = 450 ml

714205 = 5 Kg

714210 = 10 kg

Our advices contained herein regarding the use of our products provide general information only without ensuring technical properties. Any influence on fabric is subject to self-determination and therefore in the users own responsibility. Our suggestions are to the best of our knowledge and no liability can be derived thereof.